



LivaNova

Advanced Circulatory Support

User-Supplied Equipment

	TandemLife (Veno-Arterial)	TandemHeart (LA-FA)	ProtekDuo (RA-PA)	TandemLung (RA-PA)	TandemLung RD (IVC-RA)
5 Fr or 6 Fr venous and arterial sheaths	✓	✓			
Transseptal sheath (Mullins, SL1, SL2 most frequently used)		✓			
Transseptal needle (BRK, RF)		✓			
Toray Inoue wire or 0.035 Amplatz super stiff, exchange length wire		✓			
4 sterile perfusion clamps	✓	✓	✓	✓	✓
2-0 silk sutures	✓	✓	✓	✓	✓
Arterial dilators (8, 10, 12, 14 Fr)	✓	✓			
0.035 Lunderquist or Amplatz super stiff exchange length wire (ProtekDuo Insertion Kit 5100-0001)			✓	✓	
PA Catheter ≤ 90 cm and 0.035 guidewire compatible (ProtekDuo Insertion Kit 5100-0001)			✓	✓	
Tubing cutter	✓	✓			
18-gauge percutaneous entry needle and scalpel	✓	✓	✓	✓	✓

CardiacAssist, Inc. dba TandemLife® | 620 Alpha Drive, Suite 2, Pittsburgh, PA 15238

24/7 Clinical Support **800 373 1607** | Customer Service **800 373 7421** | Visit **TANDEMLIFE.COM**



[Indications for use](#) | [Safety Information](#)

IM-7100124-ACS