

# Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Balloon Tip, Plastic Tubing, Female Luer Lock Connector, 5.0 mm, CP-21005



The Coronary Artery Perfusion Cannulae with or without balloon are indicated for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

<b>ORDER CODE</b>	CP-21005
<b>SIZE</b>	5.0 mm
<b>AVAILABLE COATED</b>	no
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes
<b>NO DOP</b>	yes