

## Venous Return Cannulae, Optiflow Venous, 29 Fr, V182-29



The Venous cannulae are intended to be used to cannulate the major venous vessels during cardiopulmonary surgery for periods of up to six hours. Optiflow Cannulae are provided with a one-piece swirled and grooved tip with multiple holes. Depending on country registration, an optional phospholipid surface treatment may be available. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative

<b>ORDER CODE</b>	V182-29
<b>SIZE</b>	29 Fr
<b>AVAILABLE COATED</b>	Depending on country registration
<b>EFFECTIVE LENGTH (CM)</b>	39
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes
<b>NO DOP</b>	yes