

Venous Return Cannulae, Optiflow Venous, 29 Fr, V182-29



The Venous cannulae are intended to be used to cannulate the major venous vessels during cardiopulmonary surgery for periods of up to six hours. Optiflow Cannulae are provided with a one-piece swirled and grooved tip with multiple holes. Depending on country registration, an optional phospholipid surface treatment may be available. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative

Order Code	V182-29
Size	29 Fr
Available Coated	Depending on country registration
Effective Length (cm)	39
Qty Per Box	10
Latex-Free	yes
No DOP	yes