

# Cardioplegia Cannulae, Antegrade Aortic Root, with and without Vent Line, 4.0 mm, R502-26



The Aortic root cannulae are intended to cannulate the aortic root for cardioplegic solution delivery into the coronary arteries, and for venting of the heart during cardiopulmonary bypass surgery for periods up to six hours. The Aortic root cannulae with pressure line are intended to cannulate the aortic root for cardioplegic solution delivery into the coronary arteries, and for venting of the heart during cardiopulmonary bypass surgery for periods up to six hours. The pressure line allows pressure measurements in the aortic root. The Coronary perfusion cannulae are intended to individually cannulate the coronary ostia and to selectively deliver cardioplegic solution during cardiopulmonary bypass surgery for periods up to six hours.

Order Code	R502-26
Size	4.0 mm
Insertion Needle	12 Ga
Available Coated	Depending on country registration
Effective Length (cm)	13
Qty Per Box	10
Latex-Free	yes
No DOP	yes