

# Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Conical Tip, Stainless Steel Shaft , Female Luer Lock Connector, 11 Fr, P618-35

Placeholder  
Image

The Aortic root cannulae are intended to cannulate the aortic root for cardioplegic solution delivery into the coronary arteries, and for venting of the heart during cardiopulmonary bypass surgery for periods up to six hours. The Aortic root cannulae with pressure line are intended to cannulate the aortic root for cardioplegic solution delivery into the coronary arteries, and for venting of the heart during cardiopulmonary bypass surgery for periods up to six hours. The pressure line allows pressure measurements in the aortic root. The Coronary perfusion cannulae are intended to individually cannulate the coronary ostia and to selectively deliver cardioplegic solution during cardiopulmonary bypass surgery for periods up to six hours.

Order Code	P618-35
Size	11 Fr
Available Coated	Depending on country registration
Qty Per Box	10
Latex-Free	yes
No DOP	yes