

Femoral Cannulae, Percutaneous Cannulae Insertion Kit (Seldinger Technique), 8, 12, 14, 18, 20 Fr, HD9370202

LivaNova

Health innovation that matters

Livanova Percutaneous Cannulae Insertion Kit , for percutaneous cannulation technique, includes an 18 Gauge needle, a 0.035" guide-wire, a blade scalpel and 8,12,14, 18 and 20 Fr dilators.

ORDER CODE	HD9370202
SIZE	8, 12, 14, 18, 20 Fr
LATEX-FREE	yes
NO DOP	yes