

# Femoral Cannulae, Percutaneous Cannulae Insertion Kit (Seldinger Technique), 8, 12, 14, 18, 20 Fr, HD9370202



Livanova Percutaneous Cannulae Insertion Kit , for percutaneous cannulation technique, includes an 18 Gauge needle, a 0.035" guide-wire, a blade scalpel and 8,12,14, 18 and 20 Fr dilators.

Order Code HD9370202

Size 8, 12, 14, 18, 20 Fr

Latex-Free yes

No DOP yes