

Cardioplegia Cannulae, Retrograde, Self-Inflating PVC Balloon, 14 Fr, RCS-12114



The Retrograde Cardioplegia Cannulae are indicated for the infusion of blood or cardioplegia solution into the coronary venous system, intraoperatively. They are available with either a self- or manually-inflating balloon

Order Code	RCS-12114
------------	-----------

Size	14 Fr
------	-------

Available Coated	no
------------------	----

Effective Length (cm)	305
-----------------------	-----

Qty Per Box	5
-------------	---

Latex-Free	yes
------------	-----

Connector	Luer lock
-----------	-----------