

MICS Cannulae, Aortic Root Long Needle, 12 Ga (9 Fr), AR-17012



Livanova Aortic Root Long Needle is specifically designed for minimally invasive cardiac surgery in lateral thoracotomy, especially mitral valve procedures. It's longer than a conventional aortic root (35 cm) which allows it to be positioned in lateral thoracotomy.

ORDER CODE	AR-17012
SIZE	12 Ga (9 Fr)
INSERTION NEEDLE	14 Ga - Green
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	31,8
QTY PER BOX	10
LATEX-FREE	yes