

# MICS Cannulae, Bi-Flow Femoral Arterial, 19 Fr, 203-019

Bi-Flow is the only bidirectional arterial cannula on the market designed to reduce the risk of limb ischemia in an easy, safe and reproducible way. Unlike traditional cannulae, Bi-Flow's patented design incorporates a unique elbow and downstream perfusion channel. The cannula is easy to insert and remove and provides simultaneous systemic and distal perfusion to the lower limb.

Order Code	203-019
------------	---------

Size	19 Fr
------	-------

Available Coated	no
------------------	----

Effective Length (cm)	26,7
-----------------------	------

Qty Per Box	5
-------------	---

Latex-Free	yes
------------	-----