

MICS Cannulae, Bi-Flow Femoral Arterial, 19 Fr, 203-019



Bi-Flow is the only bidirectional arterial cannula on the market designed to reduce the risk of limb ischemia in an easy, safe and reproducible way.

Unlike traditional cannulae, Bi-Flow's patented design incorporates a unique elbow and downstream perfusion channel. The cannula is easy to insert and remove and provides simultaneous systemic and distal perfusion to the lower limb.

Order Code 203-019

Size 19 Fr

Available Coated no

Effective Length (cm) 26,7

Qty Per Box 5

Latex-Free yes