

MICS Cannulae, Bi-Flow Femoral Arterial, 19 Fr, 203-019



Bi-Flow is the only bidirectional arterial cannula on the market designed to reduce the risk of limb ischemia in an easy, safe and reproducible way.

Unlike traditional cannulae, Bi-Flow's patented design incorporates a unique elbow and downstream perfusion channel. The cannula is easy to insert and remove and provides simultaneous systemic and distal perfusion to the lower limb.

ORDER CODE	203-019
SIZE	19 Fr
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	26,7
QTY PER BOX	5
LATEX-FREE	yes