

## MICS Cannulae, Easy Flow Arterial, 23 Fr, 103-300



Livanova EasyFlow Arterial Cannulae are designed for safe and reproducible direct aortic cannulation, with a malleable stylet able to help cannula introduction through lateral thoracotomy. The configuration of the dispersion tip reduces shear stress on the aortic wall and risks of aortic plaque dislodgment. Available in the femoral configuration with a dedicated kit for Seldinger insertion.

<b>ORDER CODE</b>	103-300
<b>SIZE</b>	23 Fr
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	28,5
<b>QTY PER BOX</b>	5
<b>LATEX-FREE</b>	yes