

## MICS Cannulae, Easy Flow DUO Femoral Arterial, 19 Fr, 103-210



Livanova EasyFlow DUO Femoral Arterial Cannulae feature a dispersion tip for gentle flow in the femoral artery. The obturator allows bloodless insertion in the arterial vessel. Available in the femoral configuration with a dedicated kit for Seldinger insertion. Validated also for venous drainage support.

|                              |         |
|------------------------------|---------|
| <b>ORDER CODE</b>            | 103-210 |
| <b>SIZE</b>                  | 19 Fr   |
| <b>AVAILABLE COATED</b>      | no      |
| <b>EFFECTIVE LENGTH (CM)</b> | 28,5    |
| <b>QTY PER BOX</b>           | 5       |
| <b>LATEX-FREE</b>            | yes     |