

MICS Cannulae, Easy Flow DUO Femoral Arterial, 19 Fr, 103-210



Livanova EasyFlow DUO Femoral Arterial Cannulae feature a dispersion tip for gentle flow in the femoral artery. The obturator allows bloodless insertion in the arterial vessel. Available in the femoral configuration with a dedicated kit for Seldinger insertion. Validated also for venous drainage support.

| | |
|------------------------------|---------|
| ORDER CODE | 103-210 |
| SIZE | 19 Fr |
| AVAILABLE COATED | no |
| EFFECTIVE LENGTH (CM) | 28,5 |
| QTY PER BOX | 5 |
| LATEX-FREE | yes |