

# Suction Products, Sumps, ATRASump, 1/4", SU-29602



Livanova Suckers and Sumps are intended for use to remove excess fluid from the surgical field. They are available with flexible or malleable tubing, with weighted metal spring wound tip that terminates in a rounded metal tip. ATRASump offers a unique, a-traumatic, weighted tip with a clear cover available over the metal coil, designed with safety in mind

<b>ORDER CODE</b>	SU-29602
<b>SIZE</b>	1/4 in
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	38
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes