

Advanced Circulatory Support, Mechanical Circulatory Support, ProtekSolo Transseptal, 21 Fr, 5140- 6221



The ProtekSolo Transseptal cannula provides a powerful cannulation strategy for unloading the left ventricle. The ProtekSolo Transseptal is inserted at the femoral vein and placed in the left atrium via a transseptal procedure; oxygenated blood is returned via the femoral artery.

Order Code	5140-6221
Size	21 Fr
Available Coated	No
Effective Length (cm)	62
Qty Per Box	1
Latex-Free	yes