

Suction Products, Fluted Tip, Metal Handle, 15 Fr, S361-50



Livanova Pediatric Suckers are intended for intra-cardiac pediatric usage to ensure better visualization of the surgical field and are available with ergonomic plastic grip or rigid tubing. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative.

ORDER CODE	S361-50
SIZE	15 Fr
AVAILABLE COATED	Depending on country registration
QTY PER BOX	10
LATEX-FREE	yes
CONNECTOR	1/4 in