

## Suction Products, Fluted Tip, Plastic Handle, 15 Fr, S353-50



Livanova Pediatric Suckers are intended for intra-cardiac pediatric usage to ensure better visualization of the surgical field and are available with ergonomic plastic grip or rigid tubing. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative.

<b>ORDER CODE</b>	S353-50
<b>SIZE</b>	15 Fr
<b>AVAILABLE COATED</b>	Depending on country registration
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes
<b>CONNECTOR</b>	1/4 in