

Arterial Cannulae, PureFlex, Curved Tip, With Suture Flange, Wire-reinforced Tubing, 20 Fr/6.5 mm, RA20-2BP



The Arterial cannulae are intended to be used as perfusion cannulae to return arterial blood from the extracorporeal circuit to the patient during cardiopulmonary surgery for periods of up to six hours. The Pureflex Arterial cannulae consist of flexible tubing (tapered to a curved or straight flow tip) made of PVC with wire reinforcement

ORDER CODE	RA20-2BP
SIZE	20 Fr/6.5 mm
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	21.6
QTY PER BOX	10
LATEX-FREE	yes