

Arterial Cannulae, PureFlex, Straight Tip, With Suture Collar, Wire-reinforced Tubing, 24 Fr/8.0 mm, RA24-1SB



The arterial cannulae are intended to be used as perfusion cannulae to return arterial blood from the extracorporeal circuit to the patient during cardiopulmonary bypass surgery for periods of up to six hours.

Order Code	RA24-1SB
Size	24 Fr/8.0 mm
Available Coated	no
Effective Length (cm)	216
Qty Per Box	10
Latex-Free	yes