

Cardioplegia Cannulae, Retrograde, Manually-Inflating Silicone Balloon, 15 Fr, RCM-14315



The Retrograde Cardioplegia Cannulae are indicated for the infusion of blood or cardioplegia solution into the coronary venous system, intraoperatively. They are available with either a self- or manually-inflating balloon

| | |
|------------------------------|-----------|
| ORDER CODE | RCM-14315 |
| SIZE | 15 Fr |
| AVAILABLE COATED | no |
| EFFECTIVE LENGTH (CM) | 30,5 |
| QTY PER BOX | 10 |
| LATEX-FREE | yes |
| CONNECTOR | Luer lock |