

MICS Cannulae, Aortic Root Long Needle, 12 Ga (9 Fr), AR-17012



Livanova Aortic Root Long Needle is specifically designed for minimally invasive cardiac surgery in lateral thoracotomy, especially mitral valve procedures. It's longer than a conventional aortic root (35 cm) which allows it to be positioned in lateral thoracotomy.

| | |
|------------------------------|---------------|
| ORDER CODE | AR-17012 |
| SIZE | 12 Ga (9 Fr) |
| INSERTION NEEDLE | 14 Ga - Green |
| AVAILABLE COATED | no |
| EFFECTIVE LENGTH (CM) | 31,8 |
| QTY PER BOX | 10 |
| LATEX-FREE | yes |