

# MICS Cannulae, Easy Flow DUO Femoral Arterial, 23 Fr, 103-310



Livanova EasyFlow DUO Femoral Arterial Cannulae feature a dispersion tip for gentle flow in the femoral artery. The obturator allows bloodless insertion in the arterial vessel. Available in the femoral configuration with a dedicated kit for Seldinger insertion. Validated also for venous drainage support.

Order Code	103-310
Size	23 Fr
Available Coated	no
Effective Length (cm)	28,5
Qty Per Box	5
Latex-Free	yes