

# MICS Cannulae, RAP Femoral Venous, 22 Fr Distal/22 Fr Proximal, 200-100



Livanova RAP Femoral Venous Cannulae are designed for improved venous drainage, with a dual stage tip which drains from superior and inferior vena cava. Available in two sizes 22/22 Fr and 23/25 Fr with full-flow drainage capability. A high-quality soft introducer features a smooth transition of the cannula to improve insertion and reduce vascular damage.

<b>ORDER CODE</b>	200-100
<b>SIZE</b>	22 Fr Distal/22 Fr Proximal
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	61,5
<b>QTY PER BOX</b>	5
<b>LATEX-FREE</b>	yes