

## MICS Cannulae, Vascular Dilator Kit, 8,12,16,20, 24 Fr, 200-120



Livanova Vascular Dilator Kit, for percutaneous cannulation technique, includes an 18 Gauge needle, a 0.035 x 180 cm guide-wire, a blade scalpel and 8,12,16,20 and 24 Fr dilators.

<b>ORDER CODE</b>	200-120
<b>SIZE</b>	8,12,16,20, 24 Fr
<b>QTY PER BOX</b>	1
<b>LATEX-FREE</b>	yes