

Suction Products, Sumps, Flexible and Malleable Sucker with Fluted Tip, 21.7 Fr, SU-12402



Suction Wands and Suction Sumps are intended to be used to remove fluids from the surgical field. All of these devices are intended to be used for 6 hours or less. The devices are indicated for use in adult and pediatric patients undergoing surgical procedures requiring cardiopulmonary bypass.

Order Code SU-12402

Size 21.7 Fr

Available Coated no

Effective Length (cm) 36.8

Qty Per Box 10

Latex-Free yes