

# Venous Return Cannulae, Single Stage, Right Angle Lighthouse Tip, Wire-reinforced Tubing, 30 Fr, RV-41030



The Venous Return Cannula is indicated for use in drainage of the superior and inferior vena cava during cardiopulmonary bypass surgery. The Catheter consists of a flexible polyvinylchloride tubing wire reinforced tapered to a flow tip. The Catheters are available with a 90-degree bended distal section with a fluted distal flow tip.

|                             |          |
|-----------------------------|----------|
| <b>ORDER CODE</b>           | RV-41030 |
| <b>SIZE</b>                 | 30 Fr    |
| <b>AVAILABLE COATED</b>     | no       |
| <b>QTY PER BOX</b>          | 10       |
| <b>LATEX-FREE</b>           | yes      |
| <b>CONNECTOR ACCEPTANCE</b> | 3/8 in   |