

# Venous Return Cannulae, Single Stage, Straight Bullet Tip, Wire-reinforced Tubing, 34 Fr, V122-34



The Venous cannulae are intended to be used to cannulate the major venous vessels during cardiopulmonary surgery for periods of up to six hours. Single Stage Venous Return Cannulae are provided with straight and right angled tip, with or without reinforced tubing. Depending on country registration, an optional phospholipid surface treatment may be available. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative

<b>ORDER CODE</b>	V122-34
<b>SIZE</b>	34 Fr
<b>AVAILABLE COATED</b>	Depending on country registration
<b>EFFECTIVE LENGTH (CM)</b>	39
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes
<b>NO DOP</b>	yes
<b>CONNECTOR ACCEPTANCE</b>	3/8 in