

Vents, PVC Open Blunt Tip, with removeable Stylet, with Two Positions Latching Device, 12 Fr, VT-31212



The Vent Catheters are indicated for use during cardiopulmonary bypass surgery to vent/drain the left ventricle.

| | |
|------------------------------|----------|
| ORDER CODE | VT-31212 |
| SIZE | 12 Fr |
| AVAILABLE COATED | no |
| EFFECTIVE LENGTH (CM) | 27,2 |
| QTY PER BOX | 10 |
| LATEX-FREE | yes |