

# Venous Return Cannulae, Single Stage, Right Angle Lighthouse Tip, Wire-reinforced Tubing, 20 Fr, RV-41020



The Venous Return Cannula is indicated for use in drainage of the superior and inferior vena cava during cardiopulmonary bypass surgery. The Catheter consists of a flexible polyvinylchloride tubing wire reinforced tapered to a flow tip. The Catheters are available with a 90-degree bended distal section with a fluted distal flow tip.

<b>ORDER CODE</b>	RV-41020
<b>SIZE</b>	20 Fr
<b>AVAILABLE COATED</b>	no
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes
<b>CONNECTOR ACCEPTANCE</b>	1/4"-3/8