

# Femoral Cannulae, Percutaneous Cannulae Insertion Kit (Seldinger Technique), 8, 12, 14, 18, 20 Fr, HD9370201

**LivaNova**

Health innovation that matters

Livanova Percutaneous Cannulae Insertion Kit , for percutaneous cannulation technique, includes an 18 Gauge needle, a 0.035" guide-wire, a blade scalpel and 8,12,14, 18 and 20 Fr dilators.

<b>ORDER CODE</b>	HD9370201
<b>SIZE</b>	8, 12, 14, 18, 20 Fr
<b>LATEX-FREE</b>	yes
<b>NO DOP</b>	yes