

# Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, 3D Tip, Stainless Steel Shaft, Female Luer Lock Connector, 12 Fr, P612-40



Livanova Antegrade Cardioplegia Coronary Ostia Perfusion Cannulae are intended for cardioplegia delivery into the coronary ostia. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a “B” before their reference code. Please, contact your local sales representative.

|                         |                                   |
|-------------------------|-----------------------------------|
| <b>ORDER CODE</b>       | P612-40                           |
| <b>SIZE</b>             | 12 Fr                             |
| <b>AVAILABLE COATED</b> | Depending on country registration |
| <b>QTY PER BOX</b>      | 10                                |
| <b>LATEX-FREE</b>       | yes                               |
| <b>NO DOP</b>           | yes                               |