

Cardioplegia Cannulae, Antegrade Coronary Ostia Perfusion, Balloon Tip, Plastic Tubing, Female Luer Lock Connector, 4.0 mm, CP-21004



The Coronary Artery Perfusion Cannulae with or without balloon are indicated for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.

| | |
|-------------------------|----------|
| ORDER CODE | CP-21004 |
| SIZE | 4.0 mm |
| AVAILABLE COATED | no |
| QTY PER BOX | 10 |
| LATEX-FREE | yes |
| NO DOP | yes |