

MICS Cannulae, RAP Femoral Venous, 23 Fr Distal/25 Fr Proximal, 200-150



Livanova RAP Femoral Venous Cannulae are designed for improved venous drainage, with a dual stage tip which drains from superior and inferior vena cava. Available in two sizes 22/22 Fr and 23/25 Fr with full-flow drainage capability. A high-quality soft introducer features a smooth transition of the cannula to improve insertion and reduce vascular damage.

| | |
|------------------------------|-----------------------------|
| ORDER CODE | 200-150 |
| SIZE | 23 Fr Distal/25 Fr Proximal |
| AVAILABLE COATED | no |
| EFFECTIVE LENGTH (CM) | 65,5 |
| QTY PER BOX | 5 |
| LATEX-FREE | yes |