

MICS Cannulae, Vascular Dilator Kit, 8,12,16,20, 24 Fr, 200-120



Livanova Vascular Dilator Kit, for percutaneous cannulation technique, includes an 18 Gauge needle, a 0.035 x 180 cm guide-wire, a blade scalpel and 8,12,16,20 and 24 Fr dilators.

| | |
|--------------------|-------------------|
| ORDER CODE | 200-120 |
| SIZE | 8,12,16,20, 24 Fr |
| QTY PER BOX | 1 |
| LATEX-FREE | yes |