

# Suction Products, Suction Wands, Rigid Fluted Tip with Metal Tube and Plastic Handle, 21.7 Fr, SU- 12202



Livanova Rigid- and Flex-Tip Suction Wands are intended for intra-cardiac use. They are made of metal suction tubes with plastic handles. The Low Profile Suction Wand has a malleable low-profile metal suction tube with a plastic handle, pre-attached tubing and clamp, and terminates with a slip-on 1/4" connector. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

|                  |          |
|------------------|----------|
| Order Code       | SU-12202 |
| Size             | 21.7 Fr  |
| Available Coated | no       |
| Qty Per Box      | 10       |
| Latex-Free       | yes      |