

Suction Products, Suction Wands, Rigid Fluted Tip with Plastic Handle, 21 Fr, S351-70



Livanova Rigid- and Flex-Tip Suction Wands are intended for intra-cardiac use. They are made of metal suction tubes with plastic handles. The Low Profile Suction Wand has a malleable low-profile metal suction tube with a plastic handle, pre-attached tubing and clamp, and terminates with a slip-on 1/4" connector. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

Order Code	S351-70
Size	21 Fr
Available Coated	Depending on country registration
Qty Per Box	10
Latex-Free	yes