

## Suction Products, Sumps, Flexible and Malleable Sucker with Fluted Tip, 21.7 Fr, SU-12402



Livanova Rigid- and Flex-Tip Suction Wands are intended for intra-cardiac use. They are made of metal suction tubes with plastic handles. The Low Profile Suction Wand has a malleable low-profile metal suction tube with a plastic handle, pre-attached tubing and clamp, and terminates with a slip-on 1/4" connector. Depending on country registration, an optional phospholipid coating may be available to help improve biocompatibility. Note on ordering: coated products have a "B" before their reference code. Please, contact your local sales representative.

<b>ORDER CODE</b>	SU-12402
<b>SIZE</b>	21.7 Fr
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	36,8
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes