

Suction Products, Sumps, Flexible and Malleable Sump with Wire-wound Tip, 1/4", SU-20802



Livanova Suckers and Sumps are intended for use to remove excess fluid from the surgical field. They are available with flexible or malleable tubing, with weighted metal spring wound tip that terminates in a rounded metal tip

ORDER CODE	SU-20802
SIZE	1/4 in
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	41,9
QTY PER BOX	10
LATEX-FREE	yes