

# Vents, PVC Open Blunt Tip, with removeable Stylet, with Two Positions Latching Device, 18 Fr, VT-53218



The Vent Catheters are indicated for use during cardiopulmonary bypass surgery to vent/drain the left ventricle.

<b>ORDER CODE</b>	VT-53218
<b>SIZE</b>	18 Fr
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	36
<b>QTY PER BOX</b>	10
<b>LATEX-FREE</b>	yes
<b>NO DOP</b>	yes