

# Arterial Cannulae, Optiflow Arterial, Curved Tip, Wire- reinforced Tubing, 21 Fr/ 7.0 mm, A282-70A



The arterial cannulae are intended to be used as perfusion cannulae to return arterial blood from the extracorporeal circuit to the patient during cardiopulmonary bypass surgery for periods of up to six hours.

Order Code                      A282-70A

Size                                21 Fr/ 7.0 mm

Available Coated              Depending on country registration

Effective Length (cm) 24

Qty Per Box 10

Latex-Free yes

No DOP yes