

# Arterial Cannulae, PureFlex, Straight Tip, With Suture Collar, Wire-reinforced Tubing, 18 Fr/6.0 mm, RA18-1SB



The Arterial cannulae are intended to be used as perfusion cannulae to return arterial blood from the extracorporeal circuit to the patient during cardiopulmonary surgery for periods of up to six hours. The Pureflex Arterial cannulae consist of flexible tubing (tapered to a curved or straight flow tip) made of PVC with wire reinforcement

|                              |              |
|------------------------------|--------------|
| <b>ORDER CODE</b>            | RA18-1SB     |
| <b>SIZE</b>                  | 18 Fr/6.0 mm |
| <b>AVAILABLE COATED</b>      | no           |
| <b>EFFECTIVE LENGTH (CM)</b> | 21.6         |
| <b>QTY PER BOX</b>           | 10           |
| <b>LATEX-FREE</b>            | yes          |