

# Arterial Cannulae, PureFlex, Straight Tip, With Suture Collar, Wire-reinforced Tubing, 22 Fr/ 7.3 mm, RA22-1SB



The arterial cannulae are intended to be used as perfusion cannulae to return arterial blood from the extracorporeal circuit to the patient during cardiopulmonary bypass surgery for periods of up to six hours.

Order Code	RA22-1SB
Size	22 Fr/ 7.3 mm
Available Coated	no
Effective Length (cm)	21.6
Qty Per Box	10
Latex-Free	yes