

Cardioplegia Cannulae, Antegrade Aortic Root, with and without Vent Line, 16 Ga (5 Fr), AR-11016



The Aortic Root Cannula is designed to deliver cardioplegia solution during cardiopulmonary bypass surgery and aspirate air from the aorta at the end of the cardiopulmonary bypass procedure. The Aortic Root Cannula with vent line is designed to deliver cardioplegia solution and vent the heart during cardiopulmonary bypass surgery.

ORDER CODE	AR-11016
SIZE	16 Ga (5 Fr)
INSERTION NEEDLE	18 Ga - Pink
AVAILABLE COATED	no
EFFECTIVE LENGTH (CM)	15
QTY PER BOX	10
LATEX-FREE	yes