

## Vents, Aortic Root Aspirating Needle, 16 Ga, AX-30321



The Air Aspirator is indicated for use in delivery of cardioplegia solution during cardiopulmonary bypass surgery and aspirate air from the aorta at the end of the cardiopulmonary bypass procedure

<b>ORDER CODE</b>	AX-30321
<b>SIZE</b>	16 Ga
<b>AVAILABLE COATED</b>	no
<b>EFFECTIVE LENGTH (CM)</b>	5,1
<b>QTY PER BOX</b>	20
<b>LATEX-FREE</b>	yes