

Cardioplegia Cannulae, Flexible Vessels, VC-12010



Livanova Vessel Cannulae are intended for antegrade selective vein graft. They are made of a soft PVC tube, with a soft and flexible distal tip, designed to perfuse vessels from 3.0 to 3.6 mm in diameter. They are available in blunt tip configuration, with or without an integral one-way-check silicone valve.

ORDER CODE	VC-12010
EFFECTIVE LENGTH (CM)	5.7
QTY PER BOX	20
LATEX-FREE	yes
NO DOP	yes