

## Cardioplegia Cannulae, Flexible Vessels, VC-22010



Livanova Vessel Cannulae are intended for antegrade selective vein graft. They are made of a soft PVC tube, with a soft and flexible distal tip, designed to perfuse vessels from 3.0 to 3.6 mm in diameter. They are available in blunt tip configuration, with or without an integral one-way-check silicone valve.

|                              |          |
|------------------------------|----------|
| <b>ORDER CODE</b>            | VC-22010 |
| <b>EFFECTIVE LENGTH (CM)</b> | 5.7      |
| <b>QTY PER BOX</b>           | 20       |
| <b>LATEX-FREE</b>            | yes      |
| <b>NO DOP</b>                | yes      |