

Cardioplegia Cannulae, Flexible Vessels, VC-22010



Livanova Vessel Cannulae are intended for antegrade selective vein graft. They are made of a soft PVC tube, with a soft and flexible distal tip, designed to perfuse vessels from 3.0 to 3.6 mm in diameter. They are available in blunt tip configuration, with or without an integral one-way-check silicone valve.

| | |
|------------------------------|----------|
| ORDER CODE | VC-22010 |
| EFFECTIVE LENGTH (CM) | 5.7 |
| QTY PER BOX | 20 |
| LATEX-FREE | yes |
| NO DOP | yes |